

Quick Start Instruction

For ZIVE SP1, ZIVE MP1, & ZIVE PP1e

ZIVE SP1 single channel



ZIVE MP1 multichannel



ZIVE PP1e portable



- ① Insert installation CD and click SmartManager button to install Smart Manager program.



- ② The Smart Manager installation program will run.

- Follow the prompts through the rest of the installation process.
- If Microsoft visual C++2008 was not installed on your PC, it will be installed. Or if VC++2008 was installed on your PC already, click cancel button.



- ③ Click Labview RTE2013 button, which is needed for data analysis packages.



- ④ Install data analysis software packages.

- Click **ZMAN** button to install EIS data analysis software package
- Click **PathFinder** button to install EIS data management tool
- Click **IVMAN** button to install DC data analysis software package
- Click **SIM4U** button to install cyclic voltammetry digital simulation software package



- ⑤ Plug in the ZIVE instrument's power cord.

- Plug in the USB cable.
- Plug in the cell cable.
- Switch on the ZIVE instrument.



Plug in the potentiostat's power cord.



Plug in the USB cable.




Plug in the cell cable.




Switch on the ZIVE instrument.

- ⑥ Microsoft Windows will detect ZIVE instrument.



Run Device Manager on control panel.


Double click on Zive Lab Devices




Click "Reinstall Driver" button.



Select "Install from a list or specific location (Advanced)".



Click "Browse" button.
Select "C:\Program Files\Zive Lab\SM1e\Driver" folder and click "Next" button.



Hardware update is completed.

- ⑦ Run sm.exe in c:\program files\zivelab\sm folder.

Note:
Software manual and hardware manual are located at c:\program files\zivelab\sm\manual folder.

ZIVE SM Smart Manager

- Standard software includes;

- Basic techniques with standard functions
- | | |
|--------------------|-----------------------|
| 1. Potentiostatic | 2. Galvanostatic |
| 3. OCP measurement | 4. Potential sweep |
| 5. Current sweep | 6. Cyclic voltammetry |

- Optional software package includes;

- EIS software package(EISe)
 - Potentiostatic EIS
 - Galvanostatic EIS
 - Pseudo Galvanostatic EIS
 - OCP EIS
 - Potentiodynamic PEIS
 - Galvanodynamic GEIS
 - Potentiostatic HFR
 - Galvanostatic HFR
 - Potentiodynamic HFR
 - Galvanodynamic HFR
 - Multisine Potentiostatic EIS
 - Multisine Galvanostatic EIS
- Battery/Energy software package(BATE)
 - CC/CV charging discharging test
 - CC/CC charging discharging test
 - Discharging test
 - EVS test
 - Variable scan rate CV
 - Pstat IV curve
 - Gstat IV curve
 - SteadyState CV
- Corrosion software package(CORE)
 - Tafel(Tafel experiment)
 - Rp(Polarization resistance)
 - RpEc trend
 - PDYN(Potentiodynamic)
 - CYPOL(Cyclic polarization resistance)
 - GDYN(Galvanodynamic)
 - Reactivoion potential
 - Ecorr vs. time
 - Galvanic corrosion
 - Potentiostatic ECN
 - Glavanostatic ECN
 - ZRA mode ECN
- Electrochemical analysis software package(EASE)
 - Step techniques
 - CA(Chronoamperometry)
 - CC(Chronocoulometry)
 - CP(Chronopotentiometry)
 - Sweep techniques
 - LSV(Linear sweep voltammetry)
 - SDV(Sampled DC voltammetry)
 - Pulsed techniques
 - DVP(Differential pulse voltammetry)
 - SWV(Square wave voltammetry)
 - DPA(Differential pulse amperometry)
 - NPV(Normal pulsed voltammetry)
 - RNPV(Reverse normal pulse voltammetry)
 - DNPV(Differential normal pulse voltammetry)


Find your solution with us!



Battery



Super Capacitor



Solar Cell



Fuel Cell



Corrosion



Sensor



General
Electrochemistry

Designed by

ZIVE LAB



WonATech Co., Ltd.

7, NeungAnMal 1-gil, Seocho-gu, Seoul, 137-180, Korea

Phone: +82-2-578-6516, Fax: +82-2-576-2635, e-mail: sales@wonatech.com

website: www.wonatech.com / www.zivelab.com

