



Because the conductivity of a material is directly linked with ohmic losses, the measurement of ionic conductivity is crucially important in order to evaluate the performance of a newly synthesized material such as ion exchange membrane(IEM) and proton exchange membrane(polymer electrolyte membrane, PEM).

Today ion exchange membranes are receiving considerable attention and are successfully applied for desalination of sea and brackish water and for treating industrial effluents. And proton exchange membrane(PEM) is one of the key components for various consumer related applications for fuel cells, e.g. automobiles, back-up power, portable power etc.

For example, in PEMs, protons can transport in two directions, across the membrane and through the membrane. This results in two conductivities, in-plane conductivity and through-plane conductivity. For PEM fuel cells, through-plane conductivity measurement is more meaningful than in-plane because proton transfer occurs in the through-plane direction.

The conductivity of the membrane can be calculated based on the measured resistance by the following equation:

$$\sigma = \frac{L}{RWT}$$

where σ is the membrane conductivity(S/cm), L is the length between the electrodes, R is the measured resistance, W is the membrane width, and T is the membrane thickness.


MCJ1 Through-Plane Conductivity Test Jig

The **MCJ1** Through-plane conductivity test jig helps user to setup a 2-probe electrochemical cell consisting of 2 stainless steel probes that sandwiches the membrane to measure through-plane conductivity of membranes. The **MCJ1** is designed to hold a membrane by pulling a lever.

Normally, a number of galvanostatic alternating current(AC) electrochemical impedance spectroscopic (EIS) techniques or DC techniques are used to estimate the membrane conductivity. User can set up a perfect system with one of ZIVE series Electrochemical Workstation with **MCJ1** conductivity test jig for through-plane conductivity measurements.



MCJ1 (Through-Plane Conductivity Test Jig)



● Applications

- Polymer Electrolyte Membrane(PEM) for Fuel Cell
- Ion Exchange Membrane for Redox Flow Battery
- Water Softening/Water Purification/Water Treatment
- Desalination
- Ion Separation, etc.


● Related Product

- Membrane Conductivity Cell(MCC)
 - designed to measure ionic conductivity by simply loading a membrane into cell hardware fixture
 - 4 point probe type
 - material
 - cell body : PEEK
 - wire : platinum
 - operating temperature : to 130°C
 - fuel cell hardware available : 5, 9 and 25 cm² fuel cell test hardware fixture
 - not included, provided by WonATech

● Specifications

| | |
|-------------------------|--------------------------------|
| Measurement type | through-plane |
| Sample size | >30mm diameter |
| Sample thickness | max. 40mm |
| Sample contact material | 304 stainless steel |
| Operating temperature | up to 60°C (70°C upon request) |
| Overall dimensions | 70 x 135 x 174 mm(WxDxH) |
| Connection | 4mm banana plug |

All specifications are subject to change without notice.



MCC(Membrane Conductivity Cell)



MCC with fuel cell hardware fixture



WonATech Co., Ltd.
7 Neunganmal 1-gil, Seocho-gu,
Seoul, 06801, Korea
Tel: +82-2-578-6516 Fax: +82-2-576-2635
e-mail: sales@wonatech.com
website: www.wonatech.com

Local Distributor