

# Product Catalog

## WBCS3000 Series Battery Cell Test System



*For  
Battery/Fuel Cell  
Supercapacitor/Solar Cell  
Corrosion  
Material Testing*

## Battery Charge/Discharge System

- Battery cycle life test
- Fuel cell test
- Supercapacitor test
- Electrochemical experiment
- Potentiostat/Galvanostat circuit
- Various safety functions
- Universal graphic function
- Other application

The **WBCS series** are ideal for battery lifetime testing and its system hardware is designed for easy channel expansion and maintenance. The **WBCS3000Ls(Le)/S/M1/M2 series** consist of plug-in type modules with independent power supplies per substation, where 8 modules(8 channels) with up to 200 Watt per channel can be installed. The Dual channel module **WBCS3000D series** is for 400Watt application. The **WBCS3000H8 series** channel can be used for powers up to 800Watt per channel and **WBCS3000H12 series** channel was for max 1200Watt per channel. The **WBCS3000HP series** can be designed for powers up to 4kWatt per channel configured inside of rack. The **WBCS3000Ls32**, **WBCS3000Le32** and **WBCS3000Lx32** are comprised of 32 modules(32 channels) per substation and channels can be expanded by unit of 16 channels.

For **WBCS3000H series** and **WBCS3000HP series**, each of the channels should have its own power supply. **WBCS3000D** has power supply per 2 channels. If the cyler system consists of multiple substations or controllers, each substation or controller can be used as independent system with optional Stand Alone Kit. And extra substations can be added for channel expansion. The current range can be configured differently for each channel. These features give users a lot of flexibility in application.

Because these multichannel battery cyler systems provide independent control over each channel, user can test multiple samples simultaneously and independently. The powerful software can give various plots for users to analyze test results easily.

The **WBCS series** are designed with 4 probe type true potentiostat/galvanostat circuit and it can support most various control parameters with various cutoff conditions. Pre-defined techniques for battery test allow user to test their battery so easily even if he/she is a novice for battery testing. And it also provides powerful test schedule function for advanced users. Consequently, user can design experiments with customized schedule with various control parameters.

The **WBCS series** can also support electrochemical techniques such as corrosion test techniques, electro-analytical techniques, cyclic voltammetry, chronoamperometry, and potentiometry, etc. and this feature allows user to perform general Echem experiments.

There are 11 kinds of **WBCS model** depending on channel power or fixed specification.

### Low Power Type

- **WBCS3000Ls**  
Voltage range:  $\pm 5V$   
Max. current:  $\pm 10mA$
- **WBCS3000Le**  
Voltage range:  $\pm 5V$   
Max. current:  $\pm 100mA$
- **WBCS3000Ls32**  
Voltage range:  $\pm 5V$   
Max. current:  $\pm 10mA$   
Plug-in channel type: 32 channel per substation
- **WBCS3000Le32**  
Voltage range:  $\pm 5V$   
Max. current:  $\pm 100mA$   
Plug-in channel type: 32 channel per substation
- **WBCS3000Lx32**  
Voltage range:  $-1V$  to  $+5V$   
Max. current:  $\pm 1A$   
Plug-in channel type: 32 channel per substation

### Mid Power Type

- **WBCS3000S**  
Max. power: 50Watt  
Max. current:  $\pm 5A$
- **WBCS3000M1**  
Max. power: 100Watt  
Max. current:  $\pm 10A$
- **WBCS3000M2**  
Max. power: 200Watt  
Max. current:  $\pm 20A$

### High Power Type

- **WBCS3000D**  
Max. power: 400Watt/ch
- **WBCS3000H8**  
Max. power: 800Watt
- **WBCS3000H12**  
Max. power: 1200Watt
- **WBCS3000HP**  
Max. power:  $<4kWatt$

Slave channel type: independent housing and power supply per channel.  
It needs an 8 channel controller separately.

## Features

- **Potentiostat/Galvanostat Circuit**
  - No switching time (charging to discharging, discharging to charging)
  - Analog feedback control to keep constant voltage & current
  - Capability of electrochemical experiments by controlling voltage versus reference electrode for positive voltage polarity only. (For  $\pm$  voltage range application, choose WMPG series multichannel potentiostat/Galvanostat)
- **High Precision**
  - 16 bit(0.0015% full scale) dual ADCs for data recording and DAC per channel for control.
  - 0.02% to 0.1% full scale voltage control and reading accuracy depending on models
  - 4 current ranges (automatic/manual selection) for WBCS3000S/Ls/Le/Lx/M1/M2/D/H8 models, 3 current ranges for WBCS3000H12 models and 1 or 3 current ranges for WBCS3000HP models. - 5 current range system is also available upon request.
  - MOSFET type linear power supply circuit
  - Shield cell cable to prevent EMI noise
  - Digital coulometer is included: Every minimum hardware sampling rate, WBCS calculates capacitance value.
- **Safety**
  - Unique "Fail check" function : To protect the system and cell itself, the experiment will stop automatically when the measured value is different from control value due to battery failure or wrong cell connection, etc. e.g. Control value: 1A, Measured value: 500mA  
The experiment will stop automatically.
  - System safety parameter : If the measured value is over system specification or user defined safety limit condition, the experiment stops automatically. : User defined safety condition setting  
User can input safety level depending on chemical properties of reactants in test cell.
  - Automatic cell connection check : Before experiment, if the cell voltage value is over than stanby voltage range, program gives the warning message for the operator to check the cell connection.
  - If operator presses stop button by mistake, a confirmation message box will appear.
  - To prevent over current, a poly-switch is located in each channel. (Low and Mid power model only)
  - Watchdog function: Stop "running channel(s)" when communication failure occurs.
  - Even if a PC failure/communication failure happens, the system will work without data loss and will store the data by max. 300,000 data points per control board.
  - If the main program is down by unstable operating system, the independent server program keeps the experiment (control & data acquisition) without dead time.
  - For high power model, emergency switch is located per channel to activate the cell off when emergency happened
- **Maintenance & System Expansion**
  - Plug-in module configures channels for the system (WBCS3000S/Ls/Le/Lx/M1/M2 models).  
: Plug-in module can be easily upgraded or can modify voltage/current range. Also, if one channel is out of order, replacing the problematic channels is the easiest way to clean the problem. Optional modules are also plug-in type. .

On the other hand, if an 8 channel system consists of only one board, user cannot use all 8 channels when one of channels is out of order

- Substation add-on type :  
The channels of WBCS series can be expanded up to 128 channels per system. Users can simply add substations when they want to expand channels. One of modular system's advantages is maintenance. In case one substation is out of order in power supply unit, the other substations still work without dead time.
- A power supply is allocated for every 8 channel substation for WBCS3000S/M1/M2/L system, whereas it is allocated for every 32 channel substation for WBCS3000Ls32/Le32/Lx32. WBCS3000D system has two channels per one power supply and WBCS3000H and WBCS3000HP has its own power supply per channel.
- Automatic data backup function
- While one or more channels is working, the calibration for the other channels in rest mode is available.
- Each channel has a poly switch instead of fuse. (Low power and Mid power system only)
- Easy calibration software
- Stable TCP/IP communication
- Automatic firmware upgrade



## Options

- **Auxiliary voltage & temperature measurement**
  - User can assign any temperature input or auxiliary voltage input to single or multiple channel(s) data set by his/her demand.
  - Customized auxiliary voltage range and temperature range is available.
- **Rack mount system available**
- **Multiple control board configuration available**
- **Battery jig**
  - Coin cell jig
  - Universal jig
  - Pouch cell jig
  - Prismatic cell jig
  - Coin cell holder

## Specifications

- Circuit: potentiostat/galvanostat linear circuit
- 4 Kelvin probe type connection
- Maximum channels per PC: 128 channels
- Voltage range: unipolar or bipolar
- Current(reading & control) accuracy : 0.02~0.1% f.s. depending on model
- Sampling time varies depending on number of channels.
  - Without option
    - 1~32 channels system: 10msec
    - 33~40 channels system: 20msec
    - 41~64 channels system: 50msec
    - 65~128 channels system: 50msec
  - With option (Aux V and/or Temperature)
    - 1~16 channels system: 10msec
    - 17~40 channels system: 20msec
    - 41~64 channels system: 50msec
    - 65~128 channels system: 50msec
- Auxiliary voltage measurement range: bipolar(option)
- Temperature measurement type: K type thermocouple(option)

## Software (Smart Interface)



Multichannel control panel



Single channel control/monitor panel


- 32bit/64bit OS environment
- TCP/IP communication
- Max. 200 steps
- Max. 10 cutoff(vertex) condition
- Max. 300,000 data point memory on control board
- Single/multichannel control panel
- Various real time plots & universal axis graphs
- Data backup function
- WYSIWYG graphics
- User friendly software

### Virtual Control Panel


• BCO (Button click operation): User can do any task just by clicking the button: NO MENU SELECTION



- Easy assignment of cycle test condition file to channel with combo box selection at anytime.
- Synchronized changes of cycle test condition for selected multiple channels.
- Real time dual channel(V & I) strip chart displays for selected channel or for all running channels with time scrolling mode or whole window mode.



Multichannel real time graph




Single channel real time graph

- Status bar displays channel status.
- Various task functions: run, stop, suspend, moving step, etc.

Step No.	Test Name	Type	Initial Start Time	Range Control	Unit Step	Cut-off Condition
1	Control	Constant	00:00:00	AuxV channel	Step Time	0.000e+00
2	Control	Constant	00:00:00	AuxV channel	Step Time	Current (dV/dt) > 1000 mV x 1000
3	Control	Constant	00:00:00	AuxV channel	Step Time	0.000e+00
4	Control	Constant	00:00:00	AuxV channel	Step Time	0.000e+00

Jump step function

- Assign temperature and AuxV channel on virtual control panel



- Single channel & multichannel control/monitor panel

• VOI(Value of Interest) bar is located in upper side of 'control/monitor' window. The first value on the far left indicates the elapsed time after the experiment starts. The other parameters can be selected or be hidden at user's discretion. The selectable parameters are voltage, current, auxiliary, temperature, load, power, capacity, and energy.


Elapsed Time	Voltage(V)	Current(A)	Power(W)	Auxiliary(V)	Temperature(°C)
00:00:50	1.5m	0.000	0.000	181.3m	0.0

Voltage control

Step Time	Voltage(V)	Current(A)	Power(W)	Auxiliary(V)	Temperature(°C)
00:00:20	3.887	-50.0m	-194.4m	237.4m	0.0

Current control

### Schedule Editor




- One stop test condition creation/modification
- Multiple input parameters available
- Max. 200 test steps
- Control parameters
  - Constant voltage, LastV
  - Constant current, LastI
  - Constant power
  - Constant load
  - C-rate
  - Voltage scanning, current scanning
  - Conditioning potential
  - Conditioning current
  - Rest
  - LastVscan
  - CstepV(Staircase Voltage Sweep)
  - CstepI(Staircase Current Sweep)
  - CC/CV, CL/CV, CP/CV, C-rate/CV
  - Id, Is control
- Step flow is defined by next step, loop and cycle.
- Cut-off conditions can be set by:
  - step time, voltage, current, dV/dt, dI/dt, cycle time, loop time, capacity, -dV, Whr, Ahr, temperature, Aux voltage, dT/dt, Eoc, stepend

No	Turn Step	Type	Condi.	Value	dt	And	...
1	Next Step	Step Time	z	30.00			

Cutoff condition


- Data sampling condition by each step : time, dV/dt, dI/dt, dT/dt, dV<sup>2</sup>/dt
- And/Or logic for cut-off condition settings

## Menu Selection (Pre-defined techniques)




## Energy Test


- CC/CV (Lithium battery) test menu
- CC/CC (NiCd(NiMH) battery) test menu
- Steady state CV
- Pstat IV curve
- Gstat IV curve
- EVS (Electrochemical voltage spectroscopy) test
- GITT (Galvanostatic intermittent titration technique) test
- PITT (Potentiostatic intermittent titration technique) test




EVS test raw data



EVS graph format (dQ/dV vs. V)



GITT test



PITT test\_raw data

## Electroanalytical Techniques

- Cyclic voltammetry
- Linear sweep voltammetry
- Chrono-amperometry
- Chrono-coulometry
- Chrono-potentiometry

## Corrosion Measurement

- Tafel plot
- Potentiodynamic
- Potentiostatic
- Galvanostatic
- Cyclic polarization
- Ecorr vs. time
- Linear polarization resistance

## Simple Monitor

Ch. No.	Channel Name	Step	Value	Unit	Status
Ch. 1	Voltage	0.0000	0.0000	V	Running
Ch. 2	Current	0.0000	0.0000	A	Running
Ch. 3	Temperature	0.0000	0.0000	°C	Running
Ch. 4	Capacity	0.0000	0.0000	Ah	Running
Ch. 5	Energy	0.0000	0.0000	Wh	Running
Ch. 6	Aux V	0.0000	0.0000	V	Running
Ch. 7	Calc V	0.0000	0.0000	V	Running
Ch. 8	Temp (C)	0.0000	0.0000	°C	Running
Ch. 9	Accumulate Wh	0.0000	0.0000	Wh	Running
Ch. 10	Accumulate Ah	0.0000	0.0000	Ah	Running
Ch. 11	Accumulate Q	0.0000	0.0000	C	Running
Ch. 12	Accumulate R	0.0000	0.0000	Ω	Running
Ch. 13	Accumulate I	0.0000	0.0000	A	Running
Ch. 14	Accumulate V	0.0000	0.0000	V	Running
Ch. 15	Accumulate E	0.0000	0.0000	Wh	Running
Ch. 16	Accumulate S	0.0000	0.0000	Wh	Running
Ch. 17	Accumulate T	0.0000	0.0000	Wh	Running
Ch. 18	Accumulate U	0.0000	0.0000	Wh	Running
Ch. 19	Accumulate V	0.0000	0.0000	Wh	Running
Ch. 20	Accumulate W	0.0000	0.0000	Wh	Running
Ch. 21	Accumulate X	0.0000	0.0000	Wh	Running
Ch. 22	Accumulate Y	0.0000	0.0000	Wh	Running
Ch. 23	Accumulate Z	0.0000	0.0000	Wh	Running
Ch. 24	Accumulate AA	0.0000	0.0000	Wh	Running
Ch. 25	Accumulate AB	0.0000	0.0000	Wh	Running
Ch. 26	Accumulate AC	0.0000	0.0000	Wh	Running
Ch. 27	Accumulate AD	0.0000	0.0000	Wh	Running
Ch. 28	Accumulate AE	0.0000	0.0000	Wh	Running
Ch. 29	Accumulate AF	0.0000	0.0000	Wh	Running
Ch. 30	Accumulate AG	0.0000	0.0000	Wh	Running
Ch. 31	Accumulate AH	0.0000	0.0000	Wh	Running
Ch. 32	Accumulate AI	0.0000	0.0000	Wh	Running
Ch. 33	Accumulate AJ	0.0000	0.0000	Wh	Running
Ch. 34	Accumulate AK	0.0000	0.0000	Wh	Running
Ch. 35	Accumulate AL	0.0000	0.0000	Wh	Running
Ch. 36	Accumulate AM	0.0000	0.0000	Wh	Running
Ch. 37	Accumulate AN	0.0000	0.0000	Wh	Running
Ch. 38	Accumulate AO	0.0000	0.0000	Wh	Running
Ch. 39	Accumulate AP	0.0000	0.0000	Wh	Running
Ch. 40	Accumulate AQ	0.0000	0.0000	Wh	Running
Ch. 41	Accumulate AR	0.0000	0.0000	Wh	Running
Ch. 42	Accumulate AS	0.0000	0.0000	Wh	Running
Ch. 43	Accumulate AT	0.0000	0.0000	Wh	Running
Ch. 44	Accumulate AU	0.0000	0.0000	Wh	Running
Ch. 45	Accumulate AV	0.0000	0.0000	Wh	Running
Ch. 46	Accumulate AW	0.0000	0.0000	Wh	Running
Ch. 47	Accumulate AX	0.0000	0.0000	Wh	Running
Ch. 48	Accumulate AY	0.0000	0.0000	Wh	Running
Ch. 49	Accumulate AZ	0.0000	0.0000	Wh	Running
Ch. 50	Accumulate BA	0.0000	0.0000	Wh	Running
Ch. 51	Accumulate BB	0.0000	0.0000	Wh	Running
Ch. 52	Accumulate BC	0.0000	0.0000	Wh	Running
Ch. 53	Accumulate BD	0.0000	0.0000	Wh	Running
Ch. 54	Accumulate BE	0.0000	0.0000	Wh	Running
Ch. 55	Accumulate BF	0.0000	0.0000	Wh	Running
Ch. 56	Accumulate BG	0.0000	0.0000	Wh	Running
Ch. 57	Accumulate BH	0.0000	0.0000	Wh	Running
Ch. 58	Accumulate BI	0.0000	0.0000	Wh	Running
Ch. 59	Accumulate BJ	0.0000	0.0000	Wh	Running
Ch. 60	Accumulate BK	0.0000	0.0000	Wh	Running
Ch. 61	Accumulate BL	0.0000	0.0000	Wh	Running
Ch. 62	Accumulate BM	0.0000	0.0000	Wh	Running
Ch. 63	Accumulate BN	0.0000	0.0000	Wh	Running
Ch. 64	Accumulate BO	0.0000	0.0000	Wh	Running
Ch. 65	Accumulate BP	0.0000	0.0000	Wh	Running
Ch. 66	Accumulate BQ	0.0000	0.0000	Wh	Running
Ch. 67	Accumulate BR	0.0000	0.0000	Wh	Running
Ch. 68	Accumulate BS	0.0000	0.0000	Wh	Running
Ch. 69	Accumulate BT	0.0000	0.0000	Wh	Running
Ch. 70	Accumulate BU	0.0000	0.0000	Wh	Running
Ch. 71	Accumulate BV	0.0000	0.0000	Wh	Running
Ch. 72	Accumulate BV	0.0000	0.0000	Wh	Running
Ch. 73	Accumulate BV	0.0000	0.0000	Wh	Running
Ch. 74	Accumulate BV	0.0000	0.0000	Wh	Running
Ch. 75	Accumulate BV	0.0000	0.0000	Wh	Running
Ch. 76	Accumulate BV	0.0000	0.0000	Wh	Running
Ch. 77	Accumulate BV	0.0000	0.0000	Wh	Running
Ch. 78	Accumulate BV	0.0000	0.0000	Wh	Running
Ch. 79	Accumulate BV	0.0000	0.0000	Wh	Running
Ch. 80	Accumulate BV	0.0000	0.0000	Wh	Running

- Real time display: time, voltage, current, channel status
- Channel status color display: charging, discharging, standby, idle, calibration

## Grouping

- Classification/grouping channels by user's purpose
- Labeling each group by operator name, chemistries etc.
- Group monitor is available by this setting
- Group control




## Detailed Monitor

- Displayed test data: status, running time, step number, cycle number, step time, current range, voltage, capacity, power, energy, Aux V, Calc V, temp, cycle file name, data file name, and file size.
- Detailed monitor type selection: All channels, running channels only, and grouped channels.
- Activated character only for running channels.


Item / Ch.	1	2	3	4	5	6	7	8
Status	Run	Run	Stand	Stand	Stand	Stand	Stand	Stand
Test Time	00:18:52	00:16:12	01:41:26	01:40:41	01:40:41	01:30:59	01:40:41	01:40:41
Step / Cycle	5 / 3	1 / 1	1 / 1	4 / 2	4 / 2	4 / 2	4 / 2	6 / 2
Step Time	00:00:41	00:16:12	01:41:26	00:17:22	00:17:23	00:07:19	00:17:21	00:09:50
Current Range	Stand	Stand	Stand	100mA	100mA	100mA	100mA	100mA
Current (A)	0.0000e+0	999.7352e-6	999.7279e-6	-37.8239e-6	-61.5079e-6	0.0000e+0	-40.0532e-6	-37.6841e-6
Voltage (V)	3.48186e+0	3.74955e+0	3.67179e+0	-41.6403e-3	42.2479e-3	194.0456e-6	-41.2195e-3	-40.4164e-3
Capacity (Ah)	15.94235e+0	209.9958e+6	1.70990e+9	-150.5432e+9	-614.1407e+9	0.0000e+0	-504.5559e+9	-103.6552e+9
Power (Watt)	0.0000e+0	3.74253e+3	3.67052e+3	-1.7059e+6	2.40199e+6	0.0000e+0	2.47356e+6	1.11091e+6
Energy (Wh)	60.7992e+0	1.02072e+9	6.32097e+9	4.35164e+9	11.2327e+9	0.0000e+0	10.801e+9	3.20095e+9
Aux (V)	0.0000e+0	0.0000e+0	0.0000e+0	0.0000e+0	0.0000e+0	0.0000e+0	0.0000e+0	0.0000e+0
Calc V	0.0000e+0	0.0000e+0	0.0000e+0	0.0000e+0	0.0000e+0	0.0000e+0	0.0000e+0	0.0000e+0
Temp (C)	0.0	0.0	0.0	0.0	0.0	0.0	0.0	0.0
Schedule Name	CYTA	CYTA	CYTA	CYTA	CYTA	CYTA	CYTA	CYTA
File Name	200801201	201001302	194021003	194021004	194021005	194021006	194021007	194021008
File Size	822 Kbytes	711 Kbytes	610 Kbytes	4 Mbytes	4 Mbytes	4 Mbytes	4 Mbytes	4 Mbytes

## Graphics


- Multiple plot format
- General graph
- Cycle graph
- General graph format



# Battery Cell Test System WBCS Series




- Parameter change without reloading the data file
- Data set On/Off: Data can be visible or invisible by selecting/deselecting the data set.
- Rest step data hidden function
- Advanced graph setting




Voltage vs. |capacity| graph

## • Cycle graph format



## • Tools


- Data Editor
  - General data report
  - Cycle data report




- \* Data editing
- \* Data filtering

## • Data Conversion

- Multiple data conversion(ASCII, Excel)




- Data file splitter by cycle number





## • General Function of Graph


- Multi-parameters
- Plot overlay: max. 20 plot
- Universal graphics: any combination of X,Y1,Y2,Y3,Y4 axis parameters
- Automatic updating plot with reloading button for running channel data
- Automatic/manual scale and polarity selection for each axis



- Cross-hair pointer by mouse click/arrow key displays coordinate values on graph
- Mouse zooming
- Density, specific value display




- Copy to clipboard function to use in other application software
- Grid on/off and dot/line selection
- ASCII file conversion or Excel file conversion of graph data only




## • Calibration

- User can calibrate channel(s) while other channels are running.



## • Data backup



## Independent Data Analysis Software




The WBCS data format can be used for independent data analysis software IVMAN™ at free of charge.

- IVMAN™ software package consists of
- IVMAN software
  - IVMAN differential analysis software
  - IVMAN photo voltaic cell analysis.
  - IVMAN Tafel analysis
  - IVMAN extractor
  - IVMAN peak find module




### IVMAN DA™ Battery Test Data Analysis Software

- Battery test data analysis
- Electrochemical voltage spectroscopy (dQ/dV vs. V)
- Voltage vs. c capacity analysis (V vs. Q)
- Cycle graph (Q vs. cycle)
- Differential voltage graph(dV/dQ vs. Q)



Measured data

V vs. Q




Cycle graph

dQ/dV vs. V



### IVMAN™ Photovoltaic Cell Analysis




- Automatic analysis of parameters
- open circuit voltage, open circuit current, max. power, efficiency
- photo induced current, diode quality factor, series resistance, etc.




### IVMAN™ Main Software

- Electrochemical analysis software
- Ideal for DC corrosion data analysis and electro-analytical data analysis
- Initial guessing function on Tafel analysis
- Automatic Tafel fitting
- Polarization resistance fitting
- 3D graph
- Find peak function
- Interpolation, differentiation, integration, etc.
- Reporting function




Tafel initial fitting

Tafel analysis result




Polarization resistance fitting

Polarization analysis result



Time graph

Find peak menu




CV graph

3D graph




### IVMAN™ Tafel Analysis

- Simple Tafel calculation




### Extractor

- Extracting data by cycle number or step
- Exporting ASCII file




### Peak Find Module

- Independent peak finding software



## WBCS3000Ls(Le) Low power Type



### Application

- For low current application
- Micro battery application
- Sensor application
- Electroanalytical application

The **WBCS3000Ls(Le)** series are designed for low current applications and can be a best choice for coin cell/micro battery studies. The potential control range is specified depending on customer's specification. Up to 8 independent channels can be installed per substation. Additional channels can be added up to a maximum of 128 channels.

The **WBCS3000Ls32**, **WBCS3000Le32** and **WBCS3000Lx32** are comprised of 32 modules(32 channels) per substation and channels can be expanded by unit of 16 channels. The system has 4 current ranges. It use a local area network(LAN) for communication with a computer.

### Options

- Temperature monitoring
- Auxiliary voltage monitoring
- Coin cell holder
- Battery jig

### Specifications

Control voltage range	±5V(standard) :WBCS3000Ls(Le)32 -1V to +5V: WBCS3000Lx32
Control current range	4 ranges
LED	Run: 1ea
Input impedance	10 <sup>12</sup> Ohm
Cell connection	4 probe type, alligator clip cables
Max. channel No.	128
Voltage accuracy	±0.02% f.s.
Current accuracy	±0.02% f.s.
<b>Voltage Control/Measurement</b>	
Full scale ranges	±5V(standard) :WBCS3000Ls(Le)32 -1V to +5V: WBCS3000Lx32
Resolution (16 bits)	0.15mV(standard) @±5V
<b>Current Control/Measurement</b>	
Full scale ranges	Max. 10mA@5V(WBCS3000Ls/Ls32) Max. 100mA@5V(WBCS3000Le/Le32) Max. 1A@-1V to +5V(WBCS3000Lx32)
Resolution	16 bit(0.0015% f.s)
Communication	TCP/IP
Sampling time	Without option - 1~32 channel system : 10msec - 33~40 channel system : 20msec - 41~64 channel system : 50msec - 65~128 channel system : 50msec With option (Aux V and/or Temperature) - 1~16 channel system : 10msec - 17~40 channel system : 20msec - 41~64 channel system : 50msec - 65~128 channel system : 50msec
Size	W350.2xD328.1xH84.2mm(WBCS3000Ls/Le) W541xD454.6xH253.8mm(WBCS3000Ls32/Le32) W541xD502xH317.3mm(WBCS3000Lx32)

All specifications are subject to change without notice.

## WBCS3000S Standard Type(Mid power type)



### Application

- Battery cycle life test
- Fuel cell test
- Solar cell test
- Supercapacitor test
- Material test

The **WBCS3000S** is a research grade battery charge/discharge test system in an 8-channel substation and each channel can be used independently or simultaneously. Maximum power of each channel is 50Watt and the system can be configured with custom specification not exceeding 50Watt in power range. The **WBCS3000S** has 4 current ranges, which is suitable for various electrochemical applications. The **WBCS3000S** is designed with a local area network (LAN) for communication with a computer.

### Options

- Temperature monitoring
- Auxiliary voltage monitoring
- Coin cell holder
- Battery jig

### Specifications

Control voltage range	±5V(standard) *1
Control current range	4 ranges
LED	Run: 1ea, Mode: 2ea
Input impedance	10 <sup>12</sup> Ohm
Cell connection	4 probe type, alligator clip cables
Max. channel No.	128
Voltage accuracy	±0.02% f.s.
Current accuracy	±0.02% f.s.
<b>Voltage Control/Measurement</b>	
Full scale ranges	±5V(standard) *1
Resolution (16 bits)	0.15mV(standard) *1
<b>Current Control/Measurement</b>	
Full scale ranges	Depending on system specification Max. 5A@±5V
Resolution	16 bit(0.0015% f.s)
Communication	TCP/IP
Sampling time	Without option - 1~32 channel system : 10msec - 33~40 channel system : 20msec - 41~64 channel system : 50msec - 65~128 channel system : 50msec With option (Aux V and/or Temperature) - 1~16 channel system : 10msec - 17~40 channel system : 20msec - 41~64 channel system : 50msec - 65~128 channel system : 50msec
Size	W447xH196xD454mm

\* 1: User can specify the voltage range within <80V for difference between high and low voltage. All specifications are subject to change without notice.



## WBCS3000M1 100Watt Mid Power Type



WBCS3000M1

### Application

- Suitable for power device application
- Fuel cell test
- Battery cycle life test
- Supercapacitor test

The **WBCS3000M1** system are designed for energy device application such as battery pack, solar module, and fuel cell stack, etc. The **WBCS3000M1** are derived from the standard WBCS series battery cyler system for higher power application with 4 current ranges and maximum power of each channel is 100Watt. The system can be configured with custom specification not exceeding its maximum power. Up to 8 independent channels can be installed per substation and extra channels can be added up to a maximum of 128 channels.

### Options

- Temperature monitoring
- Battery jigs
- Auxiliary voltage monitoring

### Specifications

Control voltage range	±5V(standard) *1
Control current range	4 ranges
LED	Run: 1ea, Mode: 2ea
Input impedance	10 <sup>12</sup> Ohm
Cell connection	4 probe type, alligator clip cables
Max. channel No.	128
Voltage accuracy	±0.02% f.s.
Current accuracy	±0.05% f.s.
<b>Voltage Control/M Measurement</b>	
Full scale ranges	±5V(standard) *1
Resolution (16 bits)	0.15mV(standard) *1
<b>Current Control/M Measurement</b>	
Full scale ranges	Depending on system specification Max. 100Watt for WBCS3000M1
Resolution	16 bit(0.0015% f.s)
Communication	TCP/IP
Sampling time	Without option - 1~32 channel system : 10msec - 33~40 channel system : 20msec - 41~64 channel system : 50msec - 65~128 channel system : 50msec With option (Aux V and/or Temperature) - 1~16 channel system : 10msec - 17~40 channel system : 20msec - 41~64 channel system : 50msec - 65~128 channel system : 50msec
Size	W447.1xD498.7xH287mm

\* 1: User can specify the voltage range within <80V for difference between high and low voltage. All specifications are subject to change without notice.

## WBCS3000M2 200Watt Mid Power Type



WBCS3000M2

### Application

- Suitable for power device application
- Fuel cell test
- Battery cycle life test
- Supercapacitor test

The **WBCS3000M2** system are designed for energy device application such as battery pack, solar module, and fuel cell stack, etc. The **WBCS3000M2** are derived from the standard WBCS series battery cyler system for higher power application with 4 current ranges and maximum power of each channel is 200Watt. The system can be configured with custom specification not exceeding its maximum power. Up to 8 independent channels can be installed per substation and extra channels can be added up to a maximum of 128 channels.

### Options

- Temperature monitoring
- Battery jigs
- Auxiliary voltage monitoring

### Specifications

Control voltage range	±5V(standard) *1
Control current range	4 ranges
LED	Run: 1ea, Mode: 2ea
Input impedance	10 <sup>12</sup> Ohm
Cell connection	4 probe type, alligator clip cables
Max. channel No.	128
Voltage accuracy	±0.02% f.s.
Current accuracy	±0.05% f.s.
<b>Voltage Control/M Measurement</b>	
Full scale ranges	±5V(standard) *1
Resolution (16 bits)	0.15mV(standard) *1
<b>Current Control/M Measurement</b>	
Full scale ranges	Depending on system specification Max. 200Watt for WBCS3000M2
Resolution	16 bit(0.0015% f.s)
Communication	TCP/IP
Sampling time	Without option - 1~32 channel system : 10msec - 33~40 channel system : 20msec - 41~64 channel system : 50msec - 65~128 channel system : 50msec With option (Aux V and/or Temperature) - 1~16 channel system : 10msec - 17~40 channel system : 20msec - 41~64 channel system : 50msec - 65~128 channel system : 50msec
Size	W464.7xD625.4xH375.6mm

\* 1: User can specify the voltage range within <80V for difference between high and low voltage. All specifications are subject to change without notice.

## WBCS3000HPC 8channel controller



### General

WBCS3000HPC is 8 channel controller for high power channel such as WMPG1000D,H8,H12. This can be used for independent channel without controller.

If the system has optional input (temperature input and/or Auxiliary voltage input), These option is located at this controller.

WBCS3000HPC has power supplies for the parts in channel hardware.

Hardware specification is same in WBCS3000HPC series to control 8channel. This can be extended upto 128channels.

Channel connection between controller and channels is via dsub connector at rear side of the controller.

### Optios

- Temperature monitoring
- Auxiliary voltage monitoring



## WBCS3000D 400Watt Dual Channel Type



### Application

- Power device application
- Battery/Super capacitor test
- Solar cell test, Fuel cell test
- For multichannel application

The WBCS3000D system is designed for Mid power application (max 400Watt/ch) such as battery, solar module, fuel cell, supercapacitor etc.

The WBCS3000D has dual channel in one housing.

Each dual channel housing has its own power supply.

Customized specification is available. Up to 8 independent channels can be connected to one 8 channel controller. Additional channels can be added up to a maximum of 128 channels.

Typical models for WBCS3000D are

- $\pm 5V$  @ 26Amp WBCS3000D\_526B
- -1V to 10V @ 23Amp WBCS3000D\_1023U
- -1V to 21V @ 14Amp WBCS3000D\_2114U
- -1V to 43V @ 7Amp WBCS3000D\_437U

### Options

- Temperature monitoring
- Auxiliary voltage monitoring


### Specifications

Control current range	4 ranges
LED	Run: 1ea, Mode: 2ea, Irange: 4ea
Input impedance	$10^{12}$ Ohm for $<10V$
Cell connection	4 probe type, alligator clip cables
Max channel no.	128
Voltage accuracy	$\pm 0.05\%$ f.s.( $<10V$ )
Current accuracy	$\pm 0.05\%$ f.s.

### Current&Voltage Control/Masurement

Resolution (16 bits)	0.0015% f.s
----------------------	-------------

### Full scale ranges



Maximum current depending on voltage range

- 1) Max 26A @  $\pm 5V$
- 2) Max 23A @ -1V~+10V
- 3) Max 19A @ -1V~+12V
- 4) Max 14A @ -1V~+21V
- 5) Max 13A @ -1V~+24V
- 6) Max 10A @ -1V~+31V
- 7) Max 7A @ -1V~+43V

Communication	TCP/IP
---------------	--------

Sampling time	Without option - 1~32 channel system : 10msec - 33~40 channel system : 20msec - 41~64 channel system : 50msec - 65~128 channel system : 50msec With option (Aux V and/or Temperature) - 1~16 channel system : 10msec - 17~40 channel system : 20msec - 41~64 channel system : 50msec - 65~128 channel system : 50msec
---------------	--

Size	W447xD505xH241mm
------	------------------

All specifications are subject to change without notice.

## WBCS3000H8 800Watt Power Channel Type



### Application

- Power device application
- Battery/Super capacitor test
- Solar cell test, Fuel cell test
- For multichannel application

The **WBCS3000H8** channel is designed for power application (max 800Watt/ch) such as battery pack, solar module, fuel cell stack, super capacitor, etc. The **WBCS3000H8** has its own power supply in channel module. This requires 8 channel controller.

Customized specification is available. Up to 8 independent channels can be connected to one 8 channel controller. Additional channels can be added up to a maximum of 128 channels.

Typical models for WBCS3000H8 are

- $\pm 5V$  @ 52Amp WBCS3000H8\_552B
- -1V to 10V @ 46Amp WBCS3000H8\_1046U
- -1V to 21V @ 29Amp WBCS3000H8\_2129U
- -1V to 43V @ 15Amp WBCS3000H8\_4315U

### Options

- Temperature monitoring
- Auxiliary voltage monitoring


### Specifications

Control current range	4 ranges
LED	Run: 1ea, Mode: 2ea, Irange: 4ea
Input impedance	$10^{12}$ Ohm for <10V
Cell connection	4 probe type, alligator clip cables
Max channel no.	128
Voltage accuracy	$\pm 0.05\%$ f.s.(<10V)
Current accuracy	$\pm 0.1\%$ f.s.

### Current&Voltage Control/Measurement

Resolution (16 bits)	0.0015% f.s
----------------------	-------------

### Full scale ranges



Maximum current depending on voltage range

- 1) Max 52A @  $\pm 5V$
- 2) Max 46A @ -1V~+10V
- 3) Max 38A @ -1V~+12V
- 4) Max 29A @ -1V~+21V
- 5) Max 26A @ -1V~+24V
- 6) Max 20A @ -1V~+31V
- 7) Max 15A @ -1V~+43V

Communication	TCP/IP
---------------	--------

Sampling time	Without option - 1~32 channel system : 10msec - 33~40 channel system : 20msec - 41~64 channel system : 50msec - 65~128 channel system : 50msec With option (Aux V and/or Temperature) - 1~16 channel system : 10msec - 17~40 channel system : 20msec - 41~64 channel system : 50msec - 65~128 channel system : 50msec
---------------	--

Size	W447xD505xH241 mm
------	-------------------

All specifications are subject to change without notice.

## WBCS3000H12 1200Watt Power Channel Type



### Application

- Power device application
- Battery/Super capacitor test
- Solar cell test, Fuel cell test
- For multichannel application

The **WBCS3000H12** channel is designed for power application (max 1200Watt/ch) such as battery pack, solar module, fuel cell stack, supercapacitor etc. The **WBCS3000H12** has its own power supply in channel module. This requires 8 channel controller.

Customized specification is available. Up to 8 independent channels can be connected to one 8 channel controller. Additional channels can be added up to a maximum of 128 channels.

Typical models for WBCS3000H12 are

- -1V to 5V @ 61Amp WBCS3000H12\_561B
- -1V to 10V @ 69Amp WBCS3000H12\_1069U
- -1V to 21V @ 43Amp WBCS3000H12\_2143U
- -1V to 31V @ 30Amp WBCS3000H12\_3130U
- -1V to 43V @ 23Amp WBCS3000H12\_4323U

### Options

- Temperature monitoring
- Auxiliary voltage monitoring


### Specifications

Control current range	3 ranges
LED	Run: 1ea, Mode: 2ea, Irange: 3ea
Input impedance	$10^{12}$ Ohm for <10V
Cell connection	4 probe type, alligator clip cables
Max channel no.	128
Voltage accuracy	$\pm 0.05\%$ f.s.(<10V)
Current accuracy	$\pm 0.1\%$ f.s.

### Current&Voltage Control/Measurement

Resolution (16 bits)	0.0015% f.s
----------------------	-------------

### Full scale ranges



Maximum current depending on voltage range

- 1) Max 61A @  $\pm 5V$
- 3) Max 69A @ -1V~+10V
- 4) Max 43A @ -1V~+21V
- 5) Max 39A @ -1V~+24V
- 6) Max 30A @ -1V~+31V
- 7) Max 23A @ -1V~+43V

Communication	TCP/IP
---------------	--------

Sampling time	Without option - 1~32 channel system : 10msec - 33~40 channel system : 20msec - 41~64 channel system : 50msec - 65~128 channel system : 50msec With option (Aux V and/or Temperature) - 1~16 channel system : 10msec - 17~40 channel system : 20msec - 41~64 channel system : 50msec - 65~128 channel system : 50msec
---------------	--

Size	W464.1xD626xH285.5mm
------	----------------------

All specifications are subject to change without notice.

## WBCS3000HP High Power Type



### Application

- For power device application
- Battery/Super capacitor test
- Solar cell test, Fuel cell test
- For multichannel application
- Rack mounted model

The **WBCS3000HP** series is designed for high power application such as battery pack, solar module, fuel cell stack, electroplating, etc. The **WBCS3000HP** is derived from the standard WBCS series battery cycler system for higher power application with 1 or 3 current ranges. The system can be configured to meet customer needs. Customized specification is available. Up to 8 independent channels can be connected to one 8 channel controller. Additional channels can be added up to a maximum of 128 channels.

### Options


- Temperature monitoring
- Auxiliary voltage monitoring

### Specifications

Control current range	3 or 1 range depending on power
LED	Run: 1ea, Mode: 2ea
Input impedance	$10^{12}$ Ohm
Cell connection	4 probe type, alligator clip cables
Max. channel No.	128
Voltage accuracy	$\pm 0.1\%$ f.s.
Current accuracy	$\pm 0.1\%$ f.s.

### Current&Voltage Control/Measurement

Resolution (16 bits)	0.0015% f.s
Full scale ranges	



Maximum current depending on voltage range

- 1) Max 150A @  $\pm 5V$
- 2) Max 192A @  $-2V \sim +12V$
- 3) Max 200A @  $-2V \sim +10V$
- 4) Max 128A @  $-2V \sim +21V$
- 5) Max 85A @  $-1V \sim +24V$
- 6) Max 95A @  $-2V \sim +31V$
- 7) Max 76A @  $-1V \sim +43V$
- 8) Max 42A @  $-1V \sim +50V$
- 9) Max 51A @  $-1V \sim +68V$
- 10) Max 45A @  $-1V \sim +77V$

Communication	TCP/IP
Sampling time	Without option - 1~32 channel system : 10msec - 33~40 channel system : 20msec - 41~64 channel system : 50msec - 65~128 channel system : 50msec With option (Aux V and/or Temperature) - 1~16 channel system : 10msec - 17~40 channel system : 20msec - 41~64 channel system : 50msec - 65~128 channel system : 50msec

Size W624xD900xH1600mm

All specifications are subject to change without notice.

## System Configuration Example

### Mid power Type - Mixture System

- 8 channel system
- : 10V, 1A(3ch) + 5V, 2A(1ch) + 7V, 100mA(4ch)
- WBCS3000S substation + channel 8set

### Mid power Type & Low power Type - Mixture System

- 8 channel system
- : 5V, 1A(5ch) + 10V, 1mA(3ch)
- WBCS3000S substation + channel 5set + Modified channel 3set

### Mid power Type & High Power Type - Mixture System

- 8 channel system
- : 5V, 1A(4ch) + 10V, 25A(2ch) + 21V, 15A(2ch)
- WBCS3000S substation + channel 4set + extension card 4ea + WBCS3000D 2set

### High Power Type & Mid Power Type - Mixture System

- 8 channel system
- : 21V, 15A(4ch) + 10V, 25A(2ch) + 5V, 10A(2ch)
- 8ch controller + WBCS3000D 4set

## Jigs

### Battery Jig & Coin Cell Jig

- For cylindrical cell and/or coin cell
- 2 pin PCB type or 4 pin Lever type depending on models
- Rack type is available




### Coin Cell Holder



### Pouch Cell Jig

- It can be directly connected to the WBCS systems.
- Lever type





WonATech Co., Ltd.  
7, Neunganmal 1-gil, Seocho-gu,  
Seoul, 06801, Korea  
Phone: +82-2-578-6516  
Fax: +82-576-2635  
e-mail: sales@wonatech.com  
website: www.wonatech.com

Local Distributor



ISO 9000 & ISO 14000 Qualified