




# Battery Jigs

## Battery & Coin Cell Jig

- easy to hold cylindrical cell and/or coin cell
- wide contact point with gold coated contact area
- 4 contact point type(Kelvin probe) is available to minimize voltage drop.
- individual channel operation is available.
- rack type is available.

Contact Type	Gap <sup>1)</sup> Adjustment	Available Battery	Part No.
<b>For WBCS Series Battery Cell Test System</b>			
4 probe type	available	cylindrical cell / coin cell	UCJ <sup>*2)</sup>
4 probe type	unavailable	coin cell	CCJ <sup>*F4</sup> <sup>2)</sup>
2 probe type	unavailable	coin cell	CCJ <sup>*F2</sup> <sup>2)</sup>
<b>For ZIVE Series<sup>4)</sup>, General Battery Cell Test System &amp; Potentiostat<sup>5)</sup></b>			
4 probe type	available	cylindrical cell / coin cell	UCJ <sup>*Z</sup> <sup>2)</sup>
4 probe type	unavailable	coin cell	CCJ <sup>*FZ</sup> <sup>2)</sup>
2 probe type	unavailable	coin cell	CCJ <sup>*FZ</sup> <sup>2)</sup>
<b>For WBCS &amp; ZIVE Series, General Battery Cell Test System &amp; Potentiostat<sup>5)</sup></b>			
4 probe type	available	cylindrical cell / coin cell	UCJ <sup>1)</sup>

1) gap between pins 2) \* : number of channels 3) for single cell  
 4) for ZIVE MP2, MP5, MP5H, and MP5HC models  
 5) general battery cyler and potentiostat having banana connector



### CCJ8F2

- jig for 8 channel coin cells
- 2 probe type coin cell jig for WBCS series battery cyler
- gap between the pins cannot be adjustable.



### Rack Type

- rack type for multichannel systems
- front panel type which can be attached to a battery cyler directly is also available.



### CCJ8FZ2

- jig for 8 channel coin cells
- 2 probe type coin cell jig for ZIVE series, general battery cyler and potentiostat having banana connectors
- gap between the pins cannot be adjustable.



### UCJ1

- jig for a single cylindrical cell or coin cell
- easy to hold a coin cell by pulling the lever
- banana connector type


## Pouch Cell Jig

### PPJ Series

- pull down contact type
- 4 contact point type(Kelvin probe)

### APJ Series

- alligator clip contact type
- for WBCS series battery cycler



Contact Type	Part No.
pull-down type	PPJ* <sup>1)</sup>
alligator clip type	APJ* <sup>1)</sup>

1) \* : number of channels

## Coin Cell Holder

### For WBCS Series

- direct connect to WBCS's cell connector
- for 2016, 2025, 2032 coin cell  
- please specify coin cell size when ordering.
- Part No.  
- for standard WBCS system : CCH2  
- for low current WBCS system : CCH2L



### For ZIVE Series

- D-SUB connector type
- for ZIVE SP1, SP2, MP1, MP2, and BP2 model
- Part No.  
- for CR2032 coin cell : CCH3-20  
- for CR2450 coin cell : CCH3-24



WonATech Co., Ltd.  
7, Neunganmal 1-gil, Seocho-gu,  
Seoul, 06801, Korea  
Phone: +82-2-578-6516  
Fax: +82-576-2635  
e-mail) sales@wonatech.com  
website: www.wonatech.com

Local Distributor

