



# Coin Cell Jig


## Coin Cell Jig

- Easy to hold coin cell
- Jig for WBCS3000L/Le/Lx and S model
- Wide contact area with spring loaded cell
- Banana socket panel is available for 3rd parties cyclor or Potentiostat/Galanostat
- Individual channel operation is available.
- Optional jig cable length: 2m, 3m, 5m, 6m.
- Max height of coin cell: 6mm.
- Max current: 1.1A

### Ordering Guide

	Part No.	Description
1st	2PD	2 Pin Pull down type
	2PL	2 Pin Lever type
	<b>2PP</b>	<b>2 Pin PCB type</b>
	4PL	4 Pin Lever type
	4PK	4 Pin Knob type
2nd	<b>CCJ</b>	<b>Coin cell jig</b>
	PCJ	Pouch cell jig
	UCJM	Mid current Universal cell jig
	UCJH	High current Universal cell jig
	PRCJ	Prismatic cell jig
3rd	<b>Channel No.</b>	<b>Channel Quantity</b>
4th	<b>Z</b>	<b>Banana connector panel</b>
	H	For High temperature Normal type
5th	4P	4ch per panel
	<b>8P</b>	<b>8ch per panel</b>
	16P	16ch per panel
	20P	20ch per panel
6th	S	S type cell connector
	<b>L</b>	<b>L type cell connector</b>
	M	M type cell connector

ex)



WonATech Co., Ltd.  
7, Neunganmal 1-gil, Seocho-gu,  
Seoul, 06801, Korea  
Phone: +82-2-578-6516  
Fax: +82-576-2635  
e-mail: sales@wonatech.com  
website: www.wonatech.com

Local Distributor

