




Knob type 4 pins Universal Jig High current

Knob type 4 Pins Universal Jig High current

- Easy to hold high capacity cylindrical using knob
- Kelvin probe type 4 pin contact
- Wide contact area with spring loaded cell
- Individual channel operation is available.
- High current ring connectors with banana sensing connectors
- Max current: 50A
- Height adjustment by knob is available. Max height: 11cm
- Voltage sense connector: 4mm banana
- Current connector: cable core or ring connector



● Ordering Guide

	Part No.	Description
1st	2PD	2 Pin Pull down type
	2PL	2 Pin Lever type
	2PP	2 Pin PCB type
	4PL	4 Pin Lever type
	4PK	4 Pin Knob type
2nd	CCJ	Coin cell jig
	PCJ	Pouch cell jig
	UCJM	Mid current Universal cell jig
	UCJH	High current Universal cell jig
	PRCJ	Prismatic cell jig
3rd	Channel No.	Channel Quantity
4th	Z	Banana connector panel
	H	For High temperature Normal type
5th	4P	4ch per panel
	8P	8ch per panel
	16P	16ch per panel
	20P	20ch per panel
6th	S	S type cell connector
	L	L type cell connector
	M	M type cell connector



WonATech Co., Ltd.
7, Neunganmal 1-gil, Seocho-gu,
Seoul, 06801, Korea
Phone: +82-2-578-6516
Fax: +82-576-2635
e-mail) sales@wonatech.com
website: www.wonatech.com

Local Distributor

