


Lever type Pouch Cell Jig

Lever type Pouch Cell Jig

- Easy to hold pouch cell using clamp lever
- Kelvin probe type 4pin contact
- Jig for WBCS3000L/Le/Lx/S and M1/M2 model
- Wide contact area with spring loaded cell
- Banana socket panel is available for 3rd parties cyclers or Potentiostat/Galanostat
- Individual channel operation is available.
- Optional jig cable length: 2m, 3m, 5m, 6m.
- Open distance to insert cell lead: 9.8mm
- Adjustable cell contact width: 10mm ~ 63.6mm
- Max current: 15A (For M type cell connector)
10A (with optional banana connector)

ex)



● Ordering Guide

	Part No.	Description
1st	2PD	2 Pin Pull down type
	2PL	2 Pin Lever type
	2PP	2 Pin PCB type
	4PL	4 Pin Lever type
	4PK	4 Pin Knob type
2nd	CCJ	Coin cell jig
	PCJ	Pouch cell jig
	UCJM	Mid current Universal cell jig
	UCJH	High current Universal cell jig
	PRCJ	Prismatic cell jig
3rd	Channel No.	Channel Quantity
4th	Z	Banana connector panel
	H	For High temperature Normal type
5th	4P	4ch per panel
	8P	8ch per panel
	16P	16ch per panel
	20P	20ch per panel
6th	S	S type cell connector
	L	L type cell connector
	M	M type cell connector



4PL-PCJ8-4P-M



4PL-PCJ4-4P-M

Optional 4ch banana socket panel



WonATech Co., Ltd.
7, Neunganmal 1-gil, Seocho-gu,
Seoul, 06801, Korea
Phone: +82-2-578-6516
Fax: +82-576-2635
e-mail: sales@wonatech.com
website: www.wonatech.com

Local Distributor

