

Light Block Chamber

Features

A chamber for spectroscopic test.
This chamber protects other external light source except testing light source and works as faradaic cage.

Ordering Guide

Black Box	Part No.
Size : 603x315x330mm(WxDxH) Cable access : 40mm dia. hole, 2ea at the bottom	BB1
Size : 297x245x250mm(WxDxH) Cable access : 40mm dia. hole, 2ea, on the side	BB1S



BB1



BB1S



WonATech Co., Ltd.
7, Neunganmal 1-gil, Seocho-gu,
Seoul, 06801, Korea
Phone: +82-2-578-6516
Fax: +82-576-2635
e-mail) sales@wonatech.com
website: www.wonatech.com

Local Distributor



ISO 9000 & ISO 14000 Qualified