

Flat Specimen Holder

Features


Flat specimen holder for corrosion test



Flat specimen holder(FSH2 & FSH15) is for corrosion cell kit CCK series. You can make corrosion specimen using plate sample with this holder more conveniently. We provide two type of flat specimen holder depending on active area .

Ordering Guide

Flat Specimen Holder	Part No.
Active area : 11.28mm dia. Sample size : 15.5mm~22mm dia. Thickness : 0.3~5.8mm	FSH2
Active area : 15mm dia. Sample size : 18.5mm~25mm dia. Thickness : 0.3~5.8mm	FSH15



WonATech Co., Ltd.
7, Neunganmal 1-gil, Seocho-gu,
Seoul, 06801, Korea
Phone: +82-2-578-6516
Fax: +82-576-2635
e-mail) sales@wonatech.com
website: www.wonatech.com

Local Distributor

