




IVMAN™

Software features


- For electrochemistry data(CV, Tafel, etc.) analysis software
- Tafel analysis function
 - Initial guessing
 - Automatic fitting
 - Corrosion rate calculation
- Linear polarization automatic fitting and analysis function
- Peak search function
- Support various graph function
- Various analysis functions as simple math, interpolation, (semi) derivative, (semi) integral, etc
- Report function
- 3rd party product data analysis (needs license code)
- Free for WonATech products data format analysis (no needs license code)




Tafel fitting




Tafel graph after analysis




CV graph




Linear polarization fitting




Linear polarization graph after analysis



Find peak



Voltage/Current/Temperature vs. Time graph



3D graph

