


# Pouch Cell Jig

## PPJ Series

- pull down contact type
- 4 contact point type(Kelvin probe)



Pull-Down Type  
PPJ Series, 4ch



Pull-Down Type  
PPJ Series, 1ch

## APJ Series

- alligator clip contact type
- for WBCS series battery cyclers



Alligator Clip Type  
APJ Series, 64ch

Contact Type	Part No.
pull-down type	PPJ* <sup>1)</sup>
alligator clip type	APJ* <sup>1)</sup>

1) \* : number of channels



WonATech Co., Ltd.  
7, Neunganmal 1-gil, Seocho-gu,  
Seoul, 06801, Korea  
Phone: +82-2-578-6516  
Fax: +82-576-2635  
e-mail) sales@wonatech.com  
website: www.wonatech.com

Local Distributor

