

# Permeation Cell Kit

The permeation cell kit, PMC1, is a spinoff of flat specimen cell kit, FCK2 series, and is designed for permeation test. A membrane or a permeation foil can be placed between two glass half cells.

Two graphite plates which can be used as counter electrode and two Luggin capillary are included as standard. Membrane and reference electrode should be ordered separately.

## Features

- Cell kit for permeation experiment
- A membrane or a permeation foil to be placed between two glass half cells should be purchased separately.
- Electrodes are not included.



*Permeation Cell Kit*



*Water Jacketed Permeation Cell Kit*

## Specifications

Sample test area	
Middle	1 cm <sup>2</sup>
Both ends	5 cm <sup>2</sup>
Dimensions	
Cell vial volume	150ml x 2 ea
Size	396 x 100 x 127mm (WxDxH)
Chemical Compatibility	
Wetted materials	Pyrex®, Polycarbonate

All specifications are subject to change without notice.

## Ordering Guide

Permeation Cell Kit	PMC1
Water Jacketed Permeation Cell Kit	WPMC1



WonATech Co., Ltd.  
 7, Neunganmal 1-gil, Seocho-gu,  
 Seoul, 06801, Korea  
 Phone: +82-2-578-6516  
 Fax: +82-576-2635  
 e-mail) sales@wonatech.com  
 website: www.wonatech.com

Local Distributor



ISO 9000 & ISO 14000 Qualified