


# Platinum Electrodes

## Features

Pt plate electrode: the rod is isolated by a glass tube suitable for corrosion test  
 length : 150mm  
 glass tube outside diameter : 6mm




## Ordering Guide for Pt plate electrodes

Electrode Type	Active Area	Part No.
Pt Plate	1cm <sup>2</sup>	PFL1
	4cm <sup>2</sup>	PFL4
	5cm <sup>2</sup>	PFL5



## Features

Pt gauze electrode: Pt gauze electrode for bulk electrolysis  
 50 mesh  
 part number : PGE



Pt gauze electrode



WonATech Co., Ltd.  
 7, Neunganmal 1-gil, Seocho-gu,  
 Seoul, 06801, Korea  
 Phone: +82-2-578-6516  
 Fax: +82-576-2635  
 e-mail) sales@wonatech.com  
 website: www.wonatech.com

Local Distributor



ISO 9000 & ISO 14000 Qualified