



The universal electrode holder, UEH1, is designed to hold various sizes of electrode. The UEH1 has 4 holes to hold electrodes and three of them have a screw to adjust its hole size. The hole size is available from 1.6mm to 10mm. The material of plate is Teflon®, which has high resistance to chemicals and its white color helps user to recognize a tiny change of samples during experiments.

Specifications

Holes	
Number of holes	4
Hole size	1.6mm dia. x 1ea 6.2mm dia. x 1ea 9.6mm dia. x 1ea 10mm dia. x 1ea
Thickness	0.2mm
Rod	
Material	stainless steel
Diameter	6mm diameter
Length	220mm

All specifications are subject to change without notice.



Ordering Guide

Universal electrode holder	UEH1
----------------------------	------



WonATech Co., Ltd.
7, Neunganmal 1-gil, Seocho-gu,
Seoul, 06801, Korea
Phone: +82-2-578-6516
Fax: +82-576-2635
e-mail) sales@wonatech.com
website: www.wonatech.com

Local Distributor

