




Lever type Coin Cell Jig

Lever type Coin Cell Jig

- Easy to hold coin cell using clamp lever
- Jig for WBCS3000L/Le/Lx and S model
- Wide contact point with spring loaded cell
- Banana socket panel is available for 3rd parties cyclers or Potentiostat/Galanostat
- Individual channel operation is available.
- Optional jig cable length: 2m, 3m, 5m, 6m.



	Part No.	Description
1st	2P	2 Pin PCB type
	2PL	2 Pin Lever type
	4PL	4 Pin Lever type
2nd	CCJ	Coin cell jig
	UCJ	Universal cell jig
3rd	Channel No.	Channel Quantity
	Z	4 cell banana connector
4th	H	For High temperature
		Normal type
	4P	4ch per panel
5th	8P	8ch per panel
	16P	16ch per panel
	20P	20ch per panel
6th	S	S type cell connector
	L	L type cell connector
	M	M type cell connector



Optional banana socket panel



WonATech Co., Ltd.
 7, Neunganmal 1-gil, Seocho-gu,
 Seoul, 06801, Korea
 Phone: +82-2-578-6516
 Fax: +82-576-2635
 e-mail) sales@wonatech.com
 website: www.wonatech.com

Local Distributor

