


DDPCH


Dual Direction High Current Pouch Cell Holder

Dual Direction High Current Pouch Cell Holder


- Easy to hold pouch cell using Knob
- Kelvin probe type 4 pin contact
- Bidirectional and unidirectional measurements available
- Epoxy material with excellent chemical resistance and heat resistance
- Customized size is available (Please contact with WonATech)



Cell size - Uni direction



Cell size - Bi direction



WonATech Co., Ltd.
7, Neunganmal 1-gil, Seocho-gu,
Seoul, 06801, Korea
Phone: +82-2-578-6516
Fax: +82-576-2635
e-mail) sales@wonatech.com
website: www.wonatech.com

Local Distributor

