



# UCJ1A Single Channel Universal Jig

## Stand alone 4 Pins Universal Jig

- Easy to hold cylindrical or coin cell using clamp lever
- Kelvin probe type 4pin contact
- Wide contact area with spring loaded cell
- Banana connectors panel is for cycler or Potentiostat/Galanostat
- Max current: 10A
- Height adjustment is available. Max height: 100mm



banana socket panel & height ruler



WonATech Co., Ltd.  
7, Neunganmal 1-gil, Seocho-gu,  
Seoul, 06801, Korea  
Phone: +82-2-578-6516  
Fax: +82-576-2635  
e-mail) sales@wonatech.com  
website: www.wonatech.com

Local Distributor



ISO 9000 & ISO 14000 Qualified