



ZJIGL

EIS Cylindrical Cell Jig

4 Pins Lever type Cylindrical Cell Jig

- *Designed for EIS measurements on cylindrical cell*
- *Easy to hold cylindrical cell using clamp lever*
- *Kelvin probe type 4 pin contact*
- *Wide contact area with spring loaded cell*
- *Max current: 10A*
- *Available cylindrical cell size: 14500,18500,16340,10440*
- *Dimension: 186(W)x35(D)x54(H)mm*



WonATech Co., Ltd.
7, Neunganmal 1-gil, Seocho-gu,
Seoul, 06801, Korea
Phone: +82-2-578-6516
Fax: +82-576-2635
e-mail) sales@wonatech.com
website: www.wonatech.com

Local Distributor

