



ZJIGM EIS Cylindrical Cell Jig

4 Pins Lever type High Current Cylindrical Cell Jig

- *Designed for EIS measurements on high current cylindrical cell*
- *Easy to hold cylindrical cell using clamp lever*
- *Kelvin probe type 4 pin contact*
- *Wide contact area with spring loaded cell*
- *Max current: 50A*
- *Available cylindrical cell size: 18650,26650,32650,21700,32700*
- *Dimension: 198(W)x57(D)x60(H)mm*



WonATech Co., Ltd.
7, Neunganmal 1-gil, Seocho-gu,
Seoul, 06801, Korea
Phone: +82-2-578-6516
Fax: +82-576-2635
e-mail) sales@wonatech.com
website: www.wonatech.com

Local Distributor

